

MINIMUM REQUIREMENTS

	Client 2D	Client 3D
X-RAY SOURCES		
Processor	Intel i5 or superior	Intel i5 or superior
RAM	4 GB	8 GB
Hard Disk	> 300 GB	> 1 TB
Graphic Card	OpenGL 2.1 or DirectX 9/11	Dedicated graphic card for DeGIS 3D
Screen resolution (Width x Height)	1600 X 1024	1920 X 1080
Operating system	Windows 10 Pro 64-bit, OS X Catalina (10.15)	

RECOMMENDED REQUIREMENTS

Processor	Quadcore 2.6 GHz+	Octacore 3.0 GHz
RAM	8 GB	16 GB
Hard Disk	> 1 TB	> 1 TB
Graphic Card	Dedicated graphics card with at least 1 GB memory	Dedicated graphic card for DeGIS 3D
Screen resolution (Width x Height)	MacOS:1680x1050; Win:1280x1024	Full HD 1920 x 1080